

# lirtmats\_quick\_start

December 4, 2025

## 1 Quick Start

This python script describes how to use LiRTMaTS python package. The input data and retention time reference files used here are in <https://github.com/wanchanglin/lirtmats/tree/master/examples/data>.

### 1.1 Setup

The users need to load python package LAMP before using LiRTMaTS. It's functions used here are for loading data set and summarising the matching results. For details, see <https://github.com/wanchanglin/lamp>.

```
[1]: import sqlite3
import pandas as pd
from lamp import anno
import lirtmats.lirtmats as rtm
```

### 1.2 Data Loading

LiRTMaTS supports text files separated by comma (,) or tab (\t). The Microsoft's XLSX is also supported, using argument `sheet_name` to indicate which sheet is used for input data. The default is 0 for the first sheet.

Here we use a small example data set with `tsv` format. This data set includes peak list and intensity data matrix. LiRTMaTS requires peak list's name, m/z value and retention time. User needs to indicate the locations of feature name, m/z value, retention time and starting points of data matrix from data. Here they are 1, 2, 3 and 4, respectively.

```
[2]: cols = [1, 2, 3, 4]
data_fn = "./data/df_pos_3.tsv" # use tsv file
df = anno.read_peak(data_fn, cols, sep='\t')
df
```

```
[2]:
```

	name	mz	rt	D121	A122	\
0	M102T899	102.034153	898.850160	1.404584e+07	3.689953e+06	
1	M102T849	102.034154	849.085350	1.473961e+07		NaN
2	M105T45	105.042677	45.353942	5.520865e+05	1.813279e+05	
3	M105T54	105.054961	54.350049	6.669635e+05	4.833251e+06	
4	M105T48_1	105.074216	47.538626	6.310113e+05		NaN

...	...	...	...	...	...
1995	M299T296	299.233645	295.569540	8.125150e+04	1.020165e+05
1996	M300T43_1	299.919504	42.832066	5.042924e+04	NaN
1997	M300T62	300.119720	62.428854	NaN	3.914945e+05
1998	M300T285_2	300.124255	285.061758	NaN	4.602130e+05
1999	M300T288	300.181271	287.944377	7.880306e+05	1.738638e+06

  

	A125	A126	A127	A128	B131 \
0	3.598363e+06	1.138875e+07	4.887524e+06	2.104782e+06	7.288258e+06
1	5.934387e+06	NaN	4.607624e+06	5.969186e+06	3.367949e+06
2	2.734923e+05	2.342655e+05	6.241395e+04	1.068277e+05	1.192451e+05
3	2.137479e+06	1.552473e+06	1.753294e+06	2.301363e+06	NaN
4	5.199302e+05	4.302566e+05	5.650141e+05	3.635406e+05	1.096530e+06

  

...	...	...	...	...	...
1995	2.209362e+05	3.557402e+05	6.039153e+05	NaN	2.330915e+05
1996	NaN	2.222376e+05	3.763288e+05	2.094474e+05	1.163715e+05
1997	5.182468e+05	7.492101e+05	1.546338e+06	5.741346e+05	9.712791e+05
1998	4.559729e+05	9.718658e+05	3.864969e+05	3.877729e+05	1.315307e+06
1999	1.113482e+06	4.063701e+06	3.788191e+06	1.201084e+06	2.988076e+06

  

	...	E214	E215	E216	H234 \
0	...	3.125203e+06	3.608369e+06	NaN	4.763811e+06
1	...	1.276006e+07	1.490770e+07	2.880142e+06	4.263577e+06
2	...	3.092946e+04	1.788324e+05	1.810794e+05	3.225256e+05
3	...	1.186390e+06	3.001167e+06	2.558921e+06	NaN
4	...	7.882748e+05	NaN	9.822090e+05	4.974403e+05

  

...	...	...	...	...	...
1995	...	3.872671e+05	1.632064e+05	7.224218e+04	3.678394e+04
1996	...	4.035525e+05	2.032260e+05	2.700920e+05	NaN
1997	...	8.554399e+05	7.431820e+05	8.878200e+05	NaN
1998	...	2.418197e+06	2.917536e+06	9.108396e+05	4.583314e+05
1999	...	2.907005e+06	3.365814e+06	2.761628e+06	1.865813e+06

  

	H235	H236	H237	H238	H239 \
0	2.281365e+06	NaN	3.404450e+06	3.720441e+06	4.539032e+05
1	NaN	NaN	NaN	4.437697e+06	6.777076e+06
2	NaN	3.734778e+05	1.935349e+05	NaN	1.094705e+05
3	NaN	1.695460e+06	NaN	1.834140e+06	1.029692e+06
4	3.604541e+05	1.340656e+06	NaN	NaN	6.020203e+05

  

...	...	...	...	...	...
1995	9.526812e+04	5.785549e+04	6.183749e+05	NaN	2.915690e+04
1996	2.675647e+05	2.695188e+05	2.750383e+05	2.882957e+05	6.720465e+04
1997	3.625514e+05	4.987110e+05	1.393237e+06	5.217566e+05	NaN
1998	4.022556e+05	2.673259e+05	NaN	NaN	8.926295e+04
1999	1.956308e+06	NaN	2.918514e+06	NaN	NaN

H240

```

0          NaN
1    6.341930e+06
2    1.946732e+05
3    4.382618e+05
4    3.597655e+05
...
1995         NaN
1996    3.352428e+05
1997    1.257126e+05
1998    2.126753e+04
1999         NaN

```

[2000 rows x 40 columns]

Data frame `df` now includes only `name`, `mz`, `rt` and intensity data matrix.

### 1.3 Retention Time Matching

To perform retention time matching, users use either default retention time library or their own reference file. The reference file must have one column: `rt_lib` which is used for retention time matching with a range or tolerance in seconds. Also the column `ion_mode` should be required for indication of positive or negative mode matching. If `ion_mode` is not included in the reference file, all rows will be used for matching.

```

[3]: ion_mode = "pos"
# ref_path = "" # if empty, use default reference file for matching
ref_path = "./data/rt_lib_202509.tsv"
ref = rtm.read_rt(ref_path, ion_mode=ion_mode)
ref

```

```

[3]:
   identifier      metabolite_name  rt_lib \
0  ACMG_aqC18_POS_0001      MS5029_Isovaleraldehyde    24.6
1  ACMG_aqC18_POS_0002  L057_Dihydroxyfumaric acid hydrate    27.0
2  ACMG_aqC18_POS_0003      L061_Benzoic acid    27.0
3  ACMG_aqC18_POS_0004      L052_Spermine    28.2
4  ACMG_aqC18_POS_0005      L021_Spermidine    30.0
...
2827  ACMG_aqC18_POS_1412      LIM3312_Cholesterol    659.4
2828  ACMG_aqC18_POS_1413      L013_5alpha-Cholestan-3-one    672.6
2829  ACMG_aqC18_POS_1414  LIM3310_5alpha-Cholest-7-en-3beta-ol    675.0
2830  ACMG_aqC18_POS_1415      L0302_5alpha-Cholestanol    681.6
2831  ACMG_aqC18_POS_1416      L045_10Z-Nonadecenoic acid    723.6

   inchikey      ion_mod
0  QPUYECUOLPXSFU-UHFFFAOYSA-N  positive
1  SEKGMJVHSBBHRD-WZHZPDAFSA-M  positive
2  DMBUODUULYCPAK-UHFFFAOYSA-N  positive
3  XDSPGKDYYRNYJI-IUPFWZBJSA-N  positive

```

```

4      HELXLJCILKEWJH-NCGAPWICSA-N  positive
...
2827  ASOSVCXGWPDUGN-UHFFFAOYSA-N  negative
2828  XQCZBXHVTFVIFE-UHFFFAOYSA-N  negative
2829  WLFXSECCHULRRO-UHFFFAOYSA-N  negative
2830  YCIMNLLNPGFGHC-UHFFFAOYSA-N  negative
2831  QIGBRXMKCJVMJ-UHFFFAOYSA-N  negative

```

[2832 rows x 5 columns]

rt\_tol is a threshold for the retention time matching window. The unit is seconds and the default value is 5.

```

[4]: rt_tol = 5
     res = rtm.comp_match_rt(df, ref, rt_tol)
     res

```

```

[4]:      id      rt      identifier \
0      M105T45  45.353942  ACMG_aqC18_POS_0280
1      M105T45  45.353942  ACMG_aqC18_POS_0281
2      M105T45  45.353942  ACMG_aqC18_POS_0282
3      M105T45  45.353942  ACMG_aqC18_POS_0283
4      M105T45  45.353942  ACMG_aqC18_POS_0284
...
150065  M300T288  287.944377  ACMG_aqC18_POS_0942
150066  M300T288  287.944377  ACMG_aqC18_POS_0943
150067  M300T288  287.944377  ACMG_aqC18_POS_0944
150068  M300T288  287.944377  ACMG_aqC18_POS_0945
150069  M300T288  287.944377  ACMG_aqC18_POS_0945

```

```

      metabolite_name  rt_lib      inchikey \
0      L0309_Asymmetric dimethylarginine  40.5  ZDLDXNCMJBOYJV-YFKPBYRVSA-N
1      MS5037_Ribonic acid gamma-lactone  40.5  DAUAQNGYDSHRET-UHFFFAOYSA-N
2      L018_L-Dihydroorotic acid  40.5  KCDXJAYRVLXPFO-UHFFFAOYSA-N
3      L030_Stachydrine  40.5  ITECRQOOEQWFPE-UHFFFAOYSA-N
4      L072_Amino adipic acid  40.5  JYPHNHPXFNEZBR-UHFFFAOYSA-N
...
150065      MS5008_Ethyl crotonate  291.0  OZWKMVRBQXNZKK-UHFFFAOYSA-N
150066      MS5032_2-Phenyl-1-propanol  291.0  DKYWVDODHFEZIM-UHFFFAOYSA-N
150067      L015_Methyl indole-3-acetate  291.0  RTIXKCRFFJGDFG-UHFFFAOYSA-N
150068      L003_Cinnamic aldehyde  291.6  FNYLWVPVPXGIIP-UHFFFAOYSA-N
150069      L003_Cinnamic aldehyde  291.6  FNYLWVPVPXGIIP-UHFFFAOYSA-N

```

```

      ion_mod  rt_range
0      positive      5
1      positive      5
2      positive      5
3      positive      5

```

```

4      positive      5
...      ...      ...
150065 negative      5
150066 negative      5
150067 negative      5
150068 positive      5
150069 negative      5

```

```
[150070 rows x 8 columns]
```

## 1.4 Summarize Results

The function `comp_summ` in package `LAMP` summarises the retention time matching.

```
[5]: sr, mr = anno.comp_summ(df, res)
```

This function combines peak table with retention time matching results and returns two results in different formats. `sr` is single row results for each peak id in peak table `df`:

```
[6]: sr
```

```

[6]:      name      mz      rt  rt_range  \
0      M100T54  100.075925  53.810924    5.0
1      M1015T254  1014.985384  253.626177    5.0
2      M101T228  101.060060  228.125403    5.0
3      M102T849  102.034154  849.085350   NaN
4      M102T899  102.034153  898.850160   NaN
...      ...      ...      ...      ...
1995      M865T700  865.244172  700.365420   NaN
1996      M919T647  918.701782  646.988220    5.0
1997  M925T237_1  924.898294  236.964462    5.0
1998      M933T267  933.410460  266.976471    5.0
1999      M934T242  933.932365  242.395371    5.0

                                     identifier  \
0      ACMG_aqC18_POS_0389::ACMG_aqC18_POS_0389::ACMG...
1      ACMG_aqC18_POS_0782::ACMG_aqC18_POS_0783::ACMG...
2      ACMG_aqC18_POS_0654::ACMG_aqC18_POS_0654::ACMG...
3                                     NaN
4                                     NaN
...                                     ...
1995                                     NaN
1996  ACMG_aqC18_POS_1407::ACMG_aqC18_POS_1407::ACMG...
1997  ACMG_aqC18_POS_0690::ACMG_aqC18_POS_0691::ACMG...
1998  ACMG_aqC18_POS_0839::ACMG_aqC18_POS_0840::ACMG...
1999  ACMG_aqC18_POS_0720::ACMG_aqC18_POS_0721::ACMG...

                                     metabolite_name  \

```

0	L0488_Maleic acid::L0488_Maleic acid::L0321_L-...
1	L0481_3-Hydroxydecanedioic acid::LIM3308_Suber...
2	MS5018_Dimethyl maleate::MS5018_Dimethyl malea...
3	NaN
4	NaN
...	...
1995	NaN
1996	L005_Vitamin K1::L005_Vitamin K1::LIM3314_Phyl...
1997	L0306_Syringic acid::L0315_ortho-Hydroxyphenyl...
1998	MS5012_Diethyl malonate::MS5019_Trimethylaceti...
1999	LIM3312_Aspartame::MS5023_Ethyl levulinate::MS...

	rt_lib \
0	48.9::48.9::49.2::49.2::50.4::50.4::50.4::50.4...
1	249.00000000000003::249.00000000000003::249.00...
2	223.2::223.2::223.8::223.8::223.8::223.8::223...
3	NaN
4	NaN
...	...
1995	NaN
1996	642.6::642.6::643.2::643.2::647.1::647.1
1997	232.2::232.2::232.2::232.2::232.2::232.2::232...
1998	262.2::262.2::262.2::262.2::262.2::262.2::262...
1999	237.6::237.6::237.6::237.6::237.6::237.6::237...

	inchikey \
0	JJVNINGBHBWJH-UHFFFAOYSA-N::JJVNINGBHBWJH-UH...
1	FVWJYYTZTCVBKE-ROUWMTJPSA-N::TVZGACDUOSZQKY-UH...
2	KIWQWJKWBHZMDT-UHFFFAOYSA-N::KIWQWJKWBHZMDT-UH...
3	NaN
4	NaN
...	...
1995	NaN
1996	ZFDIRQKJPRINOQ-HYXAFXHYSA-N::ZFDIRQKJPRINOQ-HY...
1997	AFBPFSWMIHJQDM-UHFFFAOYSA-N::OISVCGZHLKNMSJ-UH...
1998	KEVYVLWNCKMXJX-UHFFFAOYSA-N::WTTJVINHCBCLGX-ZD...
1999	MBDOYVRWFFCFHM-SNAWJCMRSA-N::XPFVYQJUAUNWIW-UH...

	ion_mod
0	positive::negative::positive::negative::positi...
1	positive::positive::negative::negative::positi...
2	positive::negative::positive::positive::positi...
3	NaN
4	NaN
...	...
1995	NaN
1996	positive::negative::positive::negative::positi...

```

1997 positive::positive::positive::positive::positi...
1998 positive::positive::positive::negative::negati...
1999 positive::positive::positive::positive::negati...

```

[2000 rows x 9 columns]

mr is multiple rows format if the match more than once from the reference file:

[7]:

```
mr
```

[7]:

	name	mz	rt	identifier \
0	M100T54	100.075925	53.810924	ACMG_aqC18_POS_0389
1	M100T54	100.075925	53.810924	ACMG_aqC18_POS_0389
2	M100T54	100.075925	53.810924	ACMG_aqC18_POS_0390
3	M100T54	100.075925	53.810924	ACMG_aqC18_POS_0390
4	M100T54	100.075925	53.810924	ACMG_aqC18_POS_0391
...	...	...	...	...
150218	M934T242	933.932365	242.395371	ACMG_aqC18_POS_0775
150219	M934T242	933.932365	242.395371	ACMG_aqC18_POS_0772
150220	M934T242	933.932365	242.395371	ACMG_aqC18_POS_0773
150221	M934T242	933.932365	242.395371	ACMG_aqC18_POS_0774
150222	M934T242	933.932365	242.395371	ACMG_aqC18_POS_0775

	metabolite_name	rt_lib \
0	L0488_Maleic acid	48.9
1	L0488_Maleic acid	48.9
2	L0321_L-Theanine	49.2
3	L0321_L-Theanine	49.2
4	L0310_Dihydrothymine	50.4
...	...	...
150218	L035_2-Methoxybenzoic acid	247.2
150219	MS5015_Phenylglyoxal	247.2
150220	MS5021_Ethyl 2-methylacetoacetate	247.2
150221	L012_Homoveratrumic acid	247.2
150222	L035_2-Methoxybenzoic acid	247.2

	inchikey	ion_mod	rt_range
0	JJVNINGBHGJH-UHFFFAOYSA-N	positive	5.0
1	JJVNINGBHGJH-UHFFFAOYSA-N	negative	5.0
2	SULYEHGGXARJS-UHFFFAOYSA-N	positive	5.0
3	SULYEHGGXARJS-UHFFFAOYSA-N	negative	5.0
4	YPTJKHVBDCRKNF-UHFFFAOYSA-N	positive	5.0
...	...	...	...
150218	RFKITWRHKUYMRJ-UHFFFAOYSA-N	positive	5.0
150219	QWIZNVHXZXRPR-WSCXOGSTSA-N	negative	5.0
150220	BHTRKEVKTKCXOH-LBSADWJPSA-N	negative	5.0
150221	SEBFKMXJBCUCAI-UHFFFAOYSA-N	negative	5.0
150222	RFKITWRHKUYMRJ-UHFFFAOYSA-N	negative	5.0

[150223 rows x 9 columns]

All of results can be saved into a `sqlite3` database and use [DB Browser for SQLite](#) to view. Or save these results in other formats, such as TSV, CSV or XLSX, separately.

```
[8]: f_save = False          # here we do NOT save results
     db_out = "test.db"
     sr_out = "test_s.tsv"
     mr_out = "test_m.tsv"
     xlsx_out = "test.xlsx"

[9]: if f_save:
     # save all results into a sqlite3 database
     conn = sqlite3.connect(db_out)
     df[["name", "mz", "rt"]].to_sql("peaklist",
                                     conn,
                                     if_exists="replace",
                                     index=False)

     mr.to_sql("anno_mr", conn, if_exists="replace", index=False)
     sr.to_sql("anno_sr", conn, if_exists="replace", index=False)

     conn.commit()
     conn.close()

     # save results into text files
     sr.to_csv(sr_out, sep="\t", index=False)
     mr.to_csv(mr_out, sep="\t", index=False)

     # save results into Excel format
     with pd.ExcelWriter(xlsx_out, mode="w", engine="openpyxl") as writer:
         sr.to_excel(writer, sheet_name="single-row", index=False)
         mr.to_excel(writer, sheet_name="multiple-row", index=False)
```

It should be noted that saving of Excel file takes much longer time than text files.

## 1.5 End User Usages

LiRTMaTS provides two computation options: command line interface(CLI) and graphical user interface (GUI).

To use GUI, you need to open a terminal and type in:

```
$ lirtmats gui
```

To use CLI, open a terminal and type in command with required arguments, something like:

```
lirtmats cli \
  --input-data "./data/df_pos_3.tsv" \
  --input-sep "tab" \
  --col-idx "1, 2, 3, 4" \
```



```

--rt-path "" \
--rt-sep "tab" \
--rt-tol "5.0" \
--ion-mode "pos" \
--save-db \
--summ-type "xlsx" \

```

Execution of this command line will produce `df_pos_3_rtm.db` and `df_pos_3_rtm.xlsx` in the directory `./data/`. If the `summ-type` is `tsv` or `csv`, files `df_pos_3_rtm_s.tsv` or `df_pos_3_rtm_s.csv` and `df_pos_3_rtm_m.tsv` or `df_pos_3_rtm_m.csv` will be saved into `./data`.

For the best practice, you can create a bash script `.sh` (Linux and MacOS) or Windows script `.bat` to contain these CLI arguments. Change parameters in these files each time when processing new data set.

For example, there are `lirtmats_cli.sh` and `lirtmats_cli.bat` in <https://github.com/wanchanglin/lirtmats/tree/master/examples>.

- For Linux and MacOS terminal:

```

$ chmod +x lirtmats_cli.sh
$ ./lirtmats_cli.sh

```

- For Windows terminal:

```

$ lirtmats_cli.bat

```

Note that if users use `xlsx` files for input data and reference file when using GUI or CLI, all data must be in the first sheet. If you use `LiRTMaTS` functions in your python scripts, there are no such requirements.